

Requirements for Submission of gDNA for Exome or Custom Capture Sequencing

Sample Preparation:

1. DNA may be extracted by your method of choice, but the procedure should include or be followed by phenol:chloroform extraction and alcohol precipitation. A protocol is included below.
2. DNA quantitation: we strongly recommend Qubit since it is specific for dsDNA while Nanodrop measurements include contaminating ssDNA, RNA, oligonucleotides and nucleotides.
3. DNA amount: 5 µg
4. Minimum DNA concentration: 20 ng/µl
5. Minimum DNA volume: 10 µl
6. Recommended buffer: 10 mM Tris, 0.1 mM EDTA, pH7.8

Sample Submission:

1. Tube specifications: 1.5 – 2 ml snap cap or screw cap tubes.
2. Sample name must be clearly written on the tube.
3. Submit both by e-mail and with the shipment a completed NISC Sample Submission Form.
4. Submit an agarose gel photo or Bioanalyzer trace for each DNA sample.
5. NISC will confirm concentration using Qubit.

Phenol:Choroform Extraction Protocol for gDNA Clean Up

1. Sample should be in a volume of at least 100 µl for extraction. If necessary dilute with 10 mM Tris, 1 mM EDTA, pH 7.8.
2. Add an equal volume of fresh phenol:choroform:isoamyl alcohol (25:24:1) (Invitrogen cat. no. 15593-031).
3. Vortex vigorously for 30 sec. Centrifuge at top speed in a microcentrifuge* for 5 min.
4. Carefully transfer the upper layer to a fresh microfuge tube. Add an equal volume of chloroform.
5. Vortex vigorously for 30 sec. Centrifuge at top speed in a microcentrifuge for 5 min.
6. Carefully transfer the upper layer to a fresh microfuge tube. Add 1/10th volume of 3 M NaOAc and 2.5 volumes of ethanol (95-100%). Mix well and let set in ice for 15 min.
7. Centrifuge at top speed in a microcentrifuge for 30 min.
8. Carefully remove supernatant. Recentrifuge for 5 sec and carefully pipet off remaining liquid not disturbing the pellet. Hint: orienting the hinge of the

- tube upwards before centrifugation puts an invisible DNA pellet on the upper wall of the tube bottom so one can avoid it when pipetting.
9. Add 0.5 ml of 70 % ethanol. Centrifuge at top speed in a microcentrifuge for 5 min.
 10. Carefully remove supernatant. Recentrifuge for 5 sec and pipet off remaining liquid.
 11. Allow to air dry 5 min. DO NOT OVER DRY GENOMIC DNA – it can be very difficult to redissolve.
 12. Dissolve DNA pellet in volume of sterile 10 mM Tris, 0.1 mM EDTA, pH7.8 to bring concentration to ~100 µg/ml.

* Phase Lock Gel Tubes (5 Prime cat. no. 2302800) can be helpful in extraction. See <https://www.5prime.com/products/nucleic-acid-purification/organic-nucleic-acid-extraction/phase-lock-gel-.aspx> for details.